

- Disposable Blood Pressure Sensor
- AAMI Specifications
- Low Cost Disposable
- Dielectric Gel Barrier
- Fully Tested & Calibrated

### **DESCRIPTION**

The Model 1620 is a fully piezoresistive silicon pressure sensor for use in invasive blood pressure monitoring. The sensor is designed to be used with automated assembly equipment and can be dropped directly into a customer's disposable blood pressure housing. The sensor is designed to meet the requirements as described in the Association for the Advancement of Medical Instrumentation (AAMI) specification for Blood Pressure Transducers.

The pressure sensor consists of a pressure sensing element mounted on a ceramic substrate. Thick film resistors on the ceramic substrate are laser trimmed for compensation and calibration. A plastic cap is attached to the ceramic substrate to provide an easy method of attachment to the customer's assembly and protection for the sensing element. A dielectric gel is placed over the sensor to provide electrical and fluid isolation.

The Model 1620 pressure sensors are batch manufactured in a 10x12 element array on a ceramic substrate (120 units per substrate). The products are shipped in anti-static shipping containers. They can also be shipped on a tape and reel. Performance characteristics and packaging can be easily tailored on a special order basis to meet the requirements of specific customers.

#### **FEATURES**

- Low Cost, Medical Applications
- Small Size and Reliable Performance
- Gel Isolation for Liquids
- Operates from 5°C to 45°C
- Compatible for Automated Assembly
- 1% Accuracy for Replacements
- 5.0 uV/V/mmHg Sensitivity
- Customization for OEM Applications

## **APPLICATIONS**

- Surgical Procedures
- Intensive Care Units
- Infusion Pumps
- Kidney Dialysis Machines
- Vacuum Assisted Birth
- Intrauterine Monitoring

## **STANDARD RANGES**

Range	mmHg
-50 to 300	•



## PERFORMANCE SPECIFICATIONS

Supply Voltage: 6.0 Vdc

Imbient Temperature: 23°C (unless otherwise sp PARAMETERS	ecified) MIN	TYP	MAX	UNITS	NOTES	
Operating Pressure Range	-50	1117	300	mmHg	NOTES	
Over Pressure	125		300	psi		
Zero Pressure Offset	-20		20	mmHg		
Sensitivity	4.95	5.00	5.05	uV/V/mmHg		
Calibration	97.5	100	102.5	mmHg	1	
Linearity and Hysteresis (-30 to 100 mmHg)	-1	.00	1	mmHg	2	
Linearity and Hysteresis (100 to 200 mmHg)	-1		1	% Output	2	
Linearity and Hysteresis (200 to 300 mmHg)	-1.5		1.5	% Output	2	
Input Impedance	1200		3200	Ω		
Output Impedance	285		315	Ω		
Output Symmetry	0.95		1.05	Ratio	3	
Supply Voltage	2	6	10	Vdc or Vac rms		
Risk Current (@ 120 Vac rms, 60Hz)			2	uA		
Warm-Up Time		5		Seconds		
Frequency Response		1200		Hz		
Offset Drift			2	mmHg	4	
Thermal Span Shift	-0.1		0.1	%/°C	5	
Thermal Offset Shift	-0.3		0.3	mmHg/°C	5	
Phase Shift (@ 5KHz)			5	Degrees		
Light Sensitivity (3000 Foot Candle)		1		mmHg		
Defibrillator withstand (400 joules)	5			Discharges	6	
Sterilization (ETO)	3			Cycles	7	
Operating Temperature	10		40	°C		
Storage Temperature	-25		+70	°C		
Operating Product Life			168	Hours		
Shelf Life	3			Years		
Dielectric Breakdown		10,000		Vdc		
Humidity (External)	10-90% (non-condensing)					
Media Interface	Dielectric Gel					

#### Notes


- 1. Output of sensor with no pressure applied and a 150K Ohms resistor placed across +SUPPLY to +OUTPUT.
- 2. Best fit straight line.
- 3. Defined as common mode symmetry between any output and supply terminal.
- 4. Over an 8 hour time period after a 10 minute warm-up.
- 5. Over operating temperature range 10-40°C with respect to 23°C.
- 6. One discharge per minute performed by customer.
- 7. Sterilization performed by customer.

## **DIMENSIONS**

**Model 1620** 


## All dimensions are in inches [millimeters].

### **BODY STYLE N**




#### \* Dimensions indicated do not include flares.

## **BODY STYLE W**



#### \* Dimensions indicated do not include flares.

#### SECTION A -A



#### Notes

1. Materials Used:

Substrate: 96% Alumina

Transducer (Die): Silicon

Die Attach Adhesive: Room Temperature Vulcanizer Lid Adhesive: Medical Grade UV Curing Adhesive Conductor And Contact Pads: Platinum-Silver Alloy Wire Bonds And Bond Pads: Gold

Resistors: Ruthenium-Based Thick Film Paste


Solder Dams: Green Glass


Protective Gel Lid: Rad-Stable Polycarbonate Resin

2. Miniscus of Gel:


Max dimension below surface A = .035"[0.89]. Max dimension above surface A = .000"[0.000].

- 3. All dimensions taken at maximum draft.
- 4. All unspecified fillets and radii are.01 5" [0.38].
- 5. All draft angles 1 \*maximum.






## STANDARD PAD CONFIGURATIONS



## **CONNECTIONS**



### ORDERING INFORMATION



The information in this sheet has been carefully reviewed and is believed to be accurate; however, no responsibility is assumed for inaccuracies. Furthermore, this information does not convey to the purchaser of such devices any license under the patent rights to the manufacturer. Measurement Specialties, Inc. reserves the right to make changes without further notice to any product herein. Measurement Specialties, Inc. makes no warranty, representation or guarantee regarding the suitability of its product for any particular purpose, nor does Measurement Specialties, Inc. assume any liability arising out of the application or use of any product or circuit and specifically disclaims any and all liability, including without limitation consequential or incidental damages. Typical parameters can and do vary in different applications. All operating parameters must be validated for each customer application by customer's technical experts. Measurement Specialties, Inc. does not convey any license under its patent rights nor the rights of others.